

| Artikel | Korngröße in µm | Beschreibung | Max. Blattstärke in mm* | Beidseitig Bestell-Nr. | Einseitig links Bestell-Nr. | Einseitig rechts Bestell-Nr. | Preis |
|-----------------------|--------------------|-----------------------|-------------------------------|---------------------------|-----------------------------------|------------------------------------|-------|
| ● Polisher | 15 | extra fine | 0,10 | 1015 | 1015 L | 1015 R | 20,- |
| ● Strip | 25 | fine | 0,15 | 1025 | 1025 L | 1025 R | 20,- |
| ● Strip | 40 | medium | 0,20 | 1040 | 1040 L | 1040 R | 20,- |
| ● Strip | 60 | coarse | 0,25 | 1060 | 1060 L | 1060 R | 20,- |
| ● Strip | 90 | extra coarse | 0,30 | 1090 | 1090 L | 1090 R | 20,- |
| ● Strip <i>NEW</i> | 125 | super extra coarse | 0,40 | | 1125 L | 1125 R | 20,- |
| ● Opener | 15 | fine | 0,10 | 1115 | | | 22,- |
| ● Opener | 40 | medium | 0,20 | 1140 | | | 22,- |

*ca. Blattstärke bei Erstanwendung der Instrumente. Es wird empfohlen immer den ASR Spalt mit der Messlehre zu kontrollieren.

Alle oscident Strips sind Einzel- oder im 5er Blister verfügbar.

Unsere Preise verstehen sich in EUR zzgl. MwSt & Versand.





68 sek. Video

oscident®

oscident GmbH
www.oscident.de

Bommersheimerweg 11a
61348 Bad Homburg v.d.H. - Germany

+49 (0) 6172 2663956
info@oscident.de

| Article | Max. grain size in μm | Description | Max. blade size in mm* | Double-sided Order No. | One-sided left Order No. | One-sided right Order No. | Price |
|---|----------------------------------|--------------------|------------------------|------------------------|--------------------------|---------------------------|-------|
|  Polisher | 15 | extra fine | 0,10 | 1015 | 1015 L | 1015 R | 20,- |
|  Strip | 25 | fine | 0,15 | 1025 | 1025 L | 1025 R | 20,- |
|  Strip | 40 | medium | 0,20 | 1040 | 1040 L | 1040 R | 20,- |
|  Strip | 60 | coarse | 0,25 | 1060 | 1060 L | 1060 R | 20,- |
|  Strip | 90 | extra coarse | 0,30 | 1090 | 1090 L | 1090 R | 20,- |
|  Strip NEW | 125 | super extra coarse | 0,40 | | 1125 L | 1125 R | 20,- |
|  Opener | 15 | fine | 0,10 | 1115 | | | 22,- |
|  Opener | 40 | medium | 0,20 | 1140 | | | 22,- |

**approx. band size during first use of the instrument. It is recommended to always control the IPR space with a measurement gauge.*

All oscident Strips are available in packages of 1 pcs and 5 pcs.

The prices quoted in the offers are in EUR and do not include VAT or shipping costs.



68 sec. video

oscident®

oscident GmbH
www.oscident.de

Bommersheimerweg 11a
61348 Bad Homburg v.d.H. - Germany

+49 (0) 6172 2663956
info@oscident.de