

Artikel	Korngröße in µm	Beschreibung	Max. Blattstärke in mm*	Beidseitig Bestell-Nr.	Einseitig links Bestell-Nr.	Einseitig rechts Bestell-Nr.	Preis
● Polisher	15	extra fine	0,10	1015	1015 L	1015 R	20,-
● Strip	25	fine	0,15	1025	1025 L	1025 R	20,-
● Strip	40	medium	0,20	1040	1040 L	1040 R	20,-
● Strip	60	coarse	0,25	1060	1060 L	1060 R	20,-
● Strip	90	extra coarse	0,30	1090	1090 L	1090 R	20,-
● Strip <i>NEW</i>	125	super extra coarse	0,40		1125 L	1125 R	20,-
● Opener	15	fine	0,10	1115			22,-
● Opener	40	medium	0,20	1140			22,-

*ca. Blattstärke bei Erstanwendung der Instrumente. Es wird empfohlen immer den ASR Spalt mit der Messlehre zu kontrollieren.

Alle oscident Strips sind Einzel- oder im 5er Blister verfügbar.

Unsere Preise verstehen sich in EUR zzgl. MwSt & Versand.



68 sek. Video

oscident®

oscident GmbH
www.oscident.de

Bommersheimerweg 11a
61348 Bad Homburg v.d.H. - Germany

+49 (0) 6172 2663956
info@oscident.de

Article	Max. grain size in μm	Description	Max. blade size in mm*	Double-sided Order No.	One-sided left Order No.	One-sided right Order No.	Price
 Polisher	15	extra fine	0,10	1015	1015 L	1015 R	20,-
 Strip	25	fine	0,15	1025	1025 L	1025 R	20,-
 Strip	40	medium	0,20	1040	1040 L	1040 R	20,-
 Strip	60	coarse	0,25	1060	1060 L	1060 R	20,-
 Strip	90	extra coarse	0,30	1090	1090 L	1090 R	20,-
 Strip NEW	125	super extra coarse	0,40		1125 L	1125 R	20,-
 Opener	15	fine	0,10	1115			22,-
 Opener	40	medium	0,20	1140			22,-

**approx. band size during first use of the instrument. It is recommended to always control the IPR space with a measurement gauge.*

All oscident Strips are available in packages of 1 pcs and 5 pcs.

The prices quoted in the offers are in EUR and do not include VAT or shipping costs.



68 sec. video

oscident®

oscident GmbH
www.oscident.de

Bommersheimerweg 11a
61348 Bad Homburg v.d.H. - Germany

+49 (0) 6172 2663956
info@oscident.de